


KN95 Masks

How to wear properly



Position the mask on the face with two hands



Full strap over the ear. Adjust the tightness and fitting position of the strap to achieve the best comfort. Adjust the upper and lower edges of the mask to ensure it fits tightly on the face.



Place fingertips at the center of the nose clip. Mold the nose clip to the shape of the nose bridge by pressing inwards while moving your finger tips down both sides.



Cover the respirator with two hands. If air overflows while vigorously exhaling air or enters the mask while vigorously inhaling air, the masks needs to be adjusted by repeating steps 2-3.



WARNING

This mask helps protect against certain particles. Misuse can lead to sickness or death. Please read the instruction on the box carefully. For more information and further assistance.